



## Prime 17B

Simplicity, high performance and affordability:  
Prime 17B brings these benefits together in a  
perfect balance of power, consumption and usability.

Equipped with a 4.3" touch panel, steam generator and  
vacuum pump to maximize steam penetration and allow  
caregivers to easily sterilize without worry, it is ideal for  
dental clinics on a budget but don't want to sacrifice  
the quality of sterilization.



Model No.	Prime 17B	
Chamber Volume	17 Liters / 4.5 Gallons	
Chamber Size (Ø x Depth)	Ø 240 x 380 (mm) / Ø 9.4 x 15 (inch)	
Pressure Vessel Material	304 Stainless Steel	
Sterilization Temperature	121 °C / 134 °C	
Overall Dimension (W x H x D)	480 x 425 x 566 (mm) / 18.9 x 16.7 x 22.3 (inch)	
Net Weight	48 kg / 105.8 lbs	
Power Supply	230V, 50/60 Hz	
Current (Max.)	10.8 A	
Air Removal	Pre-Vacuum	
Drying	Post-Vacuum	
Air Filter	≤ 0.3 µm HEPA Filter	
Data Recorder	Built-in memory of the last 10 cycles of data Use SD card as internal storage (max. 8 GB)	
Data Printer	RS 232 port (optional thermal printer)	
Water Reservoir Capacity	4.2 Liters / 1.1 Gallons	
Waste Water Reservoir Capacity	2.5 Liters / 0.7 Gallons	
Water Consumption	820 ml	
Altitude	Under 1,000 m	
Program	<ul style="list-style-type: none"> <li>• 121 °C Regular</li> <li>• 134 °C Regular</li> <li>• 134 °C Dental Mode</li> <li>• 134 °C Flash</li> <li>• 134 °C Unwrapped</li> <li>• 134 °C Prion</li> <li>• Helix Test</li> <li>• Leakage Test</li> <li>• Dry Only</li> </ul>	

### Standard Accessories



#### Tray Set

Square tray holder with a capacity of three wire-mesh trays, made of 304 stainless steel.  
Tray Dimensions (W x H x D) :  
178 x 21.5 x 305 mm



#### Starter kit

- User Manual x 1 pc.,
- Quick Drain Tube x 1 pc.,
- Cham-mate x 2 pcs.

**Sturdy**  
Welllell Group

No. 168, Sec. 1, Zhongxing Road, Wugu District,  
New Taipei City 248012, Taiwan  
T: +886 2 2981 5245  
F: +886 2 2987 8208 (Office)  
E: foreign@sturdy.com.tw  
[www.sturdy.com.tw](http://www.sturdy.com.tw)



Sturdy reserves the rights to improve the equipment shown without prior notice.

