



# Super Microm Series

SA-260MA-R & SA-260MA

Asia's first FDA 510(k) cleared small tabletop autoclave with an intuitive user interface that provides easy access, effective and convenient sterilization.

The front-filling water design, is convenient for users sterilization process just filling water from front side to water tank.




Reliable Sterilization Partner



Model NO.	SA-260MA-R	SA-260MA
Chamber Volume	24 Liters / 6.3 Gallons	
Chamber Size	Ø 260 x 450 (mm)	
Diameter x Depth	Ø 10.2 x 17.7 (inch)	
Overall Dimensions	533 x 442 x 655 (mm)	
W x H x D	21 x 17.4 x 25.8 (inch)	
Water Reservoir Capacity	4.2 Liters / 1.1 Gallons	
Water Consumption	280 - 350 ml	
Water Filling	Manual the water tank automatic filling water into the chamber	
Net Weight	58 kg / 127.9 lbs	
Power Supply	120 V, 50/60 Hz	
Current	12 A	
Fuses	15 A x 2	
Control System	Microcontrollers	
Air Removal	Dynamic-air-removal	
Drying	Post-Vacuum	
Air Filter	≤ 0.3 µm HEPA Filter	
Data Recorder	Use USB disk as internal storage	NA
Altitude	Under 1,000 m	
Standards and Directives	FDA 510(k) K231489 ANSI/AAMI ST55 Table-Top Steam Sterilizers UL 61010-1 Safety Requirements for measurement control and laboratory use IEC/EN 61010-1 Safety Requirements for measurement control and laboratory use CAN/CSA-C22.2 No. 61010-1-12 IEC 61010-2-040 Safety requirements for sterilizers used to treat medical materials IEC 61326-1 for EMC Requirements for electrical software ASME Boiler & Pressure Vessel Code, Section VIII, Division 1 ISO 13485:2016 Quality Management System - Medical Devices	

## Standard Accessories



### Tray Set

Round tray rack, with a capacity of three trays and a tray holder made of 304 stainless steel.



### Heater Cover

A 304 stainless steel support surface.



### Starter kit

- User Manual x 1 pc.,
- Cham-mate x 2 pcs.



No. 168, Sec. 1, Zhongxing Road, Wugu District,  
 New Taipei City 248012, Taiwan  
 T: +886 2 2981 5245  
 F: +886 2 2987 8208 (Office)  
 E: foreign@sturdy.com.tw  
[www.sturdy.com.tw](http://www.sturdy.com.tw)



Sturdy reserves the rights to improve the equipment shown without prior notice.

