



Topster

**E.M.S. & DISPOSABLE
SERIES PRODUCTS**

2022

E.M.S. & DISPOSABLE SERIES PRODUCTS



Reliable Partner Sturdy

From the establishment in 1992, Sturdy has been dedicated to bring our clients the most advanced products. With the mission to be your reliable partner, Sturdy isn't only committed to the development of progressive technology, but also possesses all the requisite certifications. ISO 13485, CE certificates and Taiwan GMP are for quality management control approval and guarantee the safety of our products.



SILICONE MANUAL RESUSCITATOR KIT (Autoclavable)

Autoclavable by gas or steam sterilization. Three different sizes - Adult, Child, and Infant. Swivel joint (360°) between patient and face mask allows unrestricted movement, easy to disassemble, and to be cleaned.



SR-001

- Adult type silicone resuscitator × 1 pc.
- # 5 silicone mask × 1 pc.
- Oxygen reservoir bag (2,700ml, PVC) × 1 pc.
- Oxygen tubing (200cm) × 1 pc.
- Air way: 80, 90, 100mm × 1 pc. each
- Paper box × 1 pc.
- CE certified



SR-002

- Child type silicone resuscitator × 1 pc.
- # 3 silicone mask × 1 pc.
- Oxygen reservoir bag (2,700ml, PVC) × 1 pc.
- Oxygen tubing (200cm) × 1 pc.
- Air way: 60, 70, 80mm × 1 pc. each
- Paper box × 1 pc.
- CE certified



SR-003

- Infant type silicone resuscitator × 1 pc.
- # 0 silicone mask × 1 pc.
- Oxygen reservoir bag (900ml, PVC) × 1 pc.
- Oxygen tubing (200cm) × 1 pc.
- Air way: 40, 50, 60mm × 1 pc. each
- Paper box × 1 pc.

PVC MANUAL RESUSCITATOR KIT (Disposable)



SR-011

- Adult type PVC resuscitator × 1 pc.
- # 5 PVC mask × 1 pc.
- Oxygen reservoir bag (1,700ml, PE) × 1 pc.
- Oxygen tubing (200cm) × 1 pc.
- Air way: 80, 90, 100mm × 1 pc. each
- Paper box × 1 pc.



SR-012

- Child type PVC resuscitator × 1 pc.
- # 3 PVC mask × 1 pc.
- Oxygen reservoir bag (1,700ml, PE) × 1 pc.
- Oxygen tubing (200cm) × 1 pc.
- Air way: 60, 70, 80mm × 1 pc. each
- Paper box × 1 pc.



SR-013

- Infant type PVC resuscitator × 1 pc.
- # 0 silicone mask × 1 pc.
- Oxygen reservoir bag (1,700ml, PE) × 1 pc.
- Oxygen tubing (200cm) × 1 pc.
- Air way: 40, 50, 60mm × 1 pc. each
- Paper box × 1 pc.

Carry Box (Optional)

Perfect for carry resuscitator kit

SR-B1: For adult size

SR-B2: For child or infant size






PVC MANUAL RESUSCITATOR KIT (Disposable) - NEW




SR-011(-E)

- Adult type PVC resuscitator × 1 pc.
- # 5 Air Cushion Mask × 1 pc. (SM-065)
- Oxygen reservoir bag (1,700ml, PVC) × 1 pc.
- Oxygen tubing (200cm) × 1 pc.
- Air way: 80, 90, 100mm × 1 pc. each
- Paper box × 1 pc.
- CE certified



SR-012(-E)

- Child type PVC resuscitator × 1 pc.
- # 3 Air Cushion Mask × 1 pc. (SM-063)
- Oxygen reservoir bag (1,700ml, PVC) × 1 pc.
- Oxygen tubing (200cm) × 1 pc.
- Air way: 60, 70, 80mm × 1 pc. each
- Paper box × 1 pc.
- CE certified



SR-013(-E)

- Infant type PVC resuscitator × 1 pc.
- # 1 Air Cushion Mask × 1 pc. (SM-061)
- Oxygen reservoir bag (1,700ml, PVC) × 1 pc.
- Oxygen tubing (200cm) × 1 pc.
- Air way: 40, 50, 60mm × 1 pc. each
- Paper box × 1 pc.
- CE certified

POCKET TYPE RESUSCITATION FACE MASK (CPR KIT)


Great protection through a one-way valve and 3M filter. Air cushion design provides good seal and comfortable fit. It can be used for adult, child, and infant.



SP-007

- Carry box × 1 pc.
- One-way valve × 1 pc.
- Air filter × 1 pc.
- Air cushion Mask × 1 pc.
- Connector for manual resuscitator × 1 pc.
- Elastic strap × 1 pc.
- CE certified

SILICONE MASK (Autoclavable)



SM-030~SM-035

- Suitable for use in anesthesia, respiratory care, and resuscitation.
- Applicable to gas and steam autoclave sterilization.
- Improvement in patients' safety and comfort.
- Size: #5, #4, #3, #2, #1, #0. (For adult, child, baby, and infant use.)
- CE certified

PVC MASK (Disposable)



SP-033~SP-035

- Suitable for use in anesthesia, respiratory care, and resuscitation.
- Size: #5, #4, #3. (For adult and child use.)



AIR CUSHION MASK



- Transparent mask with cushion design provides better fitting and observation.
- Suitable for use in anesthesia, respiratory care, and resuscitation.
- Optional adjustable valve for air inflation.
- Color-coded hook ring for simple size identification.
- Size: #5, #4, #3, #2, #1, #0.
(For adult, child, baby, and infant use.)
- CE certified



- SM-060~SM-065: With valve
- SM-060A~SM-065A: No valve
- SM-260~SM-265: With vertical valve



VENTURI KIT - FIXED TYPE



SVM-1

- Mask with metal nose clip, elastic strap, 200cm tubing
- Humidity cup and reservoir
- Six interchangeable color-coded Venturi valves for different concentrations: 24%, 28%, 31%, 35%, 40%, and 50%
- Available with child size mask

VENTURI KIT - CONVERTIBLE TYPE




SVM-2

- A high-flow oxygen therapy device can supply both low and medium oxygen concentrations from 24% ~ 50%
- Mask with metal nose clip, elastic strap, 200cm tubing
- Multi-concentration valve with seven scales: 24%, 26%, 28%, 31%, 35%, 40%, and 50%.
- Available with child size mask



OXYGEN MASK – MEDIUM CONCENTRATION



Soft, transparent, and odor-free vinyl, with metal nose clip and elastic strap for patient comfort and best fit. (Made by PVC plastic)

SM-051(-E): Medium concentration mask with 210cm tubing for adult (CE certified)

SM-052(-E): Medium concentration mask with 210cm tubing for child (CE certified)

SM-053(-E): Nasal oxygen cannula (over ear) with 210cm tubing for adult (CE certified)

Two types of nasal tips at choice. (Child-sized nasal tips are also available.)

■ Standard type:

SM-053(-E) / SM-053-C(-E) (For child)

■ Curved type (Injection molding):

SM-053A(-E) / SM-053A-C(-E) (For child)

OXYGEN MASK – HIGH CONCENTRATION (NON RE-BREATHING TYPE)



Soft, transparent, and odor-free vinyl, with metal nose clip and elastic strap for patient comfort and best fit. High concentration masks with one-way valves, reservoir bag, and 210cm tubing for adult and child. (Made by PVC plastic)

SM-H051(-E): 3 one-way valves - Two valves on the mask to keep high concentration oxygen in the mask; one valve between mask and reservoir to stop exhaled gas penetrating reservoir bag. (CE certified)

SM-H051-C(-E) (For child)

NEBULIZER KITS (Disposable)



SNB-030A



SNB-030B



SNB-030C

A small volume nebulizer suitable for any position. Transparent design with scale for easy observation.

SNB-030A : Nebulizer with T-type adapter, mouthpiece, reservoir, and oxygen tube

SNB-030B : Nebulizer with oxygen mask (Adult) and oxygen tube

SNB-030C : Nebulizer with oxygen mask (Pediatric) and oxygen tube

LARGE VOLUME NEBULIZER TRACHEOTOMY MASK (Disposable)



SNB-031



STM-001



STM-002

Disposable large volume nebulizer bottle w/ adjustable oxygen concentration device, heater adapter, and DISS connector. The bottle comes with scale for max. and min. volume.

Soft, transparent, odor free vinyl, with elastic strap and swivel adapter (360°) for patient comfort and best fit. (Tubing adapter is also available)

STM-001: Adult type mask with 22mm O.D. connector

STM-002: Child type mask with 22mm O.D. connector



JAW OPENER



SM-054

Breathing relief for a trismus jaw. T-shaped for easy handling and greater leverage.

HUMIDIFIER (Disposable)



SNB-032

Unique diffuser designed to offer patients with comfortable humidity and minimal noise during therapy. The bottle is 350ml capacity, comes with scale for max. and min. volume, and 3 PSI audible alarm


MANUAL EMERGENCY SUCTION UNIT (Disposable)



SSM-01

Our compact and portable manual operating suction pump is perfect for emergency use as it is easy to carry and operate. The extra collection container and tube are included.

RIGID GUEDEL AIRWAY



Made of PE material, color-coded design with a channel in the center. (CE certified)

Model No.	Dimension (mm)	Color
SG-004(-E)	40	● Pink
SG-005(-E)	50	● Blue
SG-006(-E)	60	● Black
SG-007(-E)	70	○ White
SG-008(-E)	80	● Green
SG-009(-E)	90	● Yellow
SG-010(-E)	100	● Red
SG-011(-E)	110	● Orange

SOFT GUEDEL AIRWAY



Made of transparent soft PVC, color-coded design with a channel in the center. (CE certified)


Model No.	Dimension (mm)	Color
SSG-004	40	● Pink
SSG-005	50	● Blue
SSG-006	60	● Black
SSG-007	70	○ White
SSG-008	80	● Green
SSG-009	90	● Yellow
SSG-010	100	● Red
SSG-012	120	● Purple



LARYNGOSCOPE KIT



Fiber optic light produces a cold brighter illumination. Blades adaptable to any GREEN SPECIFICATION fiber optic system. High tensile strength handle and disposable blades suitable for C.P.R. and anesthetic operation procedure. Packed with individual pouch.




OPTIC FIBER BLADE - NEW TYPE

It is made by a whole new mold, fiber optic is fully covered to improve light efficiency and can fit THE ENTIRE green system handle much easier.


Feature :

- High impact material.
- High transparency bright light tubing, and do not generate heat.
- Ergonomic shape design for easy intubation.
- Smooth blade edge design for better patient comfort and safety.
- Light weight and compact.



Model No.	Product Description	A (mm)	B (mm)	C (mm)
MIL#0 (1110500)	Miller Style Blade	21.7	79.9	21.7
MIL#1 (1110501)	Miller Style Blade	21.7	103.9	21.7
MIL#2 (1110512)	Miller Style Blade	19.5	153.9	27.5
MIL#3 (1110513)	Miller Style Blade	19.5	194.2	28.5
MAC#1 (1110511)	Macintosh Style Blade	21	95.5	35
MAC#2 (1110502)	Macintosh Style Blade	24.9	113.9	35.7
MAC#3 (1110503)	Macintosh Style Blade	25	133.9	38.1
MAC#4 (1110504)	Macintosh Style Blade	25	152.6	38.5





Sturdy

Wellell Group

Sturdy Industrial Co., Ltd.

No.168, Sec. 1, Zhongxing Road, Wugu District,
New Taipei City 248012, Taiwan

T +886 2 2981 5245

F +886 2 2987 8208

foreign@sturdy.com.tw

www.sturdy.com.tw



Distributed by :

