



E.N.T.
Equipment



**2021
CATALOGUE**



Sturdy ENT Treatment tables offer diversity in functionality, size, and accessories to match and fit most clinical ENT applications.

SN-801/SN-802/SN-803 Specialty Treatment Tables



Features :

- All tables are Stainless Steel #304.
- The spray gun pump is oil-less and maintenance-free.
- The suction system is automatic control, and the pump is oil-less and maintenance-free.
- Throat mirror anti-fog device, timer: 10 seconds (SN-805 excluded).
- All tables fit with flexible free-moving light arms.
- Eustachian tube catheter equips with a flow control switch.
- All models include convenient free-moving medicine vial shelves.



SN-801 - 4 spray guns

SN-803 - 5 spray guns
(Artificial Marble Countertop)



Artificial Marble Countertop



SN-802 - 5 spray guns

SN-805

Small Treatment Table



SN-805
For small areas
3 spray guns

Specifications


Model	SN-801	SN-802	SN-803	SN-805
Over All Dimensions (L x W x H)	100 x 50 x 86 cm	100 x 50 x 86 cm	130 x 60 x 84 cm	76 x 46 x 85 cm
Housing Stainless Steel #304	Yes	Yes	Yes [Artificial marble countertop]	Yes
Spray Gun	4 PCs	5 PCs	5 pcs	3 PCs
Spray Gun Pump (Oil-less, Maint. free)	Yes	Yes	Yes	Yes
Suction Pump (Oil-less, Maint. free)	Yes	Yes	Yes	Yes
Moveable Free Arm Light	Yes	Yes	Yes	Yes
Eustachian Tube Catheter Device	Electrical (Optional)	Electrical (Optional)	Electrical (Optional)	N/A
Power Source	110V, 220V, or 240V 50/60 Hz			
Power Consumption	110V / 10A, 220V / 5A, or 240V / 5A			110V / 6A, or 220V / 3A, or 240V / 3A
Anti-Fog Device (350W)	Yes	Yes	Yes	N/A

Accessories

Medicine Bottle (standard)	12 PCs	12 PCs	12 PCs	8 PCs
Moveable Instrument Plate (standard)	1 PC	1 PC	1 PC	1 PC
Moveable Instrument Cup (standard)	2 PCs (12cm)	2 PCs (12cm)	2 PCs (12cm)	1 PC (10cm)
Cotton Case (standard)	4 PC	4 PC	4 PCs	2 PCs

SN-601/SN-602

Electric treatment chair



SN-601
Manual
backrest adjustment

SN-602
Electric hydraulic pump
backrest adjustment

Features :

- Multiple ways to control the treatment chair :
 - **A** Control lifting and inclination by button (SN-602)
 - **B** Control lifting by a foot pedal (SN-601 & SN-602)
 - **C** Control inclination by hand crank (SN-601)
 - **D** Control rotation by hand crank (SN-602)
- The backrest, headrest, and seat are all adjustable.

Specifications

Model	SN-602	SN-601
Dimensions (L x W x H)	76 x 62 x 56 - 84 cm	75 x 62 x 52 - 80 cm
Backrest inclination	0° ~ 80°	
Backrest Adjustment	Hydraulic Pump	Manual
Seat Rotation	90° Left and Right	N/A
Mattress	PVC foam covered with PVC plastic	
Lifting	56 ~ 84 cm	52 ~ 80 cm
Lifting Adjustment	Hydraulic Pump	
Headrest Adjustment	Manual	
Power Source	110V, 220V, or 240V 50/60Hz	
Power Consumption	110V / 7A, 220V / 4A, or 240V / 4A	



STURDY INDUSTRIAL CO., LTD.

No.168, Sec. 1, Zhongxing Road, Wugu District,
New Taipei City 24872, Taiwan

T +886 - 2 2981 - 5245

F +886 - 2 2984 - 8208

foreign@sturdy.com.tw

www.sturdy.com.tw



Sturdy reserves the right to improve the products herein without prior notice.