

Hydraulic Universal Head-End Controlled Operating Table Model No.

ST-08

Recommend You The Economical & Reliable Operating Table



STURDY INDUSTRIAL CO., LTD.

MORE SAFETY - RELIABILITY - SATISFACTION

STURDY
APEX GROUP

HYDRAULIC UNIVERSAL OPERATING TABLE ST-08

Description

STURDY designed ST-08 to meet different surgical requirements. It's a universal table at a cost-effective price.

General Design

- * Stable base construction with four castors for easy movement
- * Manual base brake
- * Lateral tilt and trendelenburg adjusted by hand crank
- * Base cover, side rails, accessories & clamps made of #304 stainless steel
- * Table top elevation controlled by hydraulic system
- * Table top was divided into 5 sections including head plate (U-shaped head rest optional), back plate, seat plate and leg plates
- * 360° table top rotation
- * With side rails for the insertion of X-ray cassettes (must equipped with acrylic sheet plates for X-ray radiography)



SPECIFICATION:

# LENGTH OF TABLE TOP	190 CM
# WIDTH OF TABLE TOP	50 CM
# TABLE TOP ELEVATION	71 - 100 CM
# LATERAL TILT LEFT/RIGHT	25° / 25°
# TRENDLENBURG / REVERSE	30° / 22°
# BACK SECTION UP/DOWN	70° / 0°
# HEAD SECTION UP/DOWN	30° / 90°
# LEG SECTION UP/DOWN	20° / 90°
# BEARING CAPACITY:	150 KGS

STANDARD ACCESSORIES:

1. SCREEN FRAME (TAC-01)	1 SET
2. SHOULDER SUPPORTER (TAC-10)	1 PR.
3. BODY SUPPORT (TAC-11)	1 PR.
4. ARM BOARD (TAC-14)	1 PR.
5. LEG CRUTCHES (TAC-13)	1 PR.
6. KIDNEY ELEVATOR (TAC-30)	1 SET

OPTIONAL ACCESSORIES:

RESTRAIN STRAP (TAC-02) / PC.
HAND GRIPS (TAC-08) / PR.
U-SHAPED HEADREST (TAC-03) / SET
RAISED ARMREST (TAC-04) / SET
X-RAY RADIOGRAPH TABLE TOP (TAC-12) / SET
CONTAMINATION BUCKET (TAC-09) / SET

* **STURDY** reserves the right to improve the equipment shown without prior notice.

STURDY
APEX GROUP

<http://www.sturdy.com.tw>

Table-ST08-2018 V.1