

Gynecology-Delivery & operating SG-680



SG-680 is a multi-purpose gynecological & obstetric operating table with electro-hydraulic controlled system. Auxiliary table storable under the table (manually articulated).

Features:

- Multifunction for sitting and lying position delivery, operation, and examination.
- Water proof foot switch and 24V control
- Full table body and accessories are covered with Stainless Steel #304, easy to keep the table clean.
- The auxiliary table is concealed under the table top, easy to extend and elevate.
- Re-moveable delivery collection plate, the table top is available for vagina doing ultrasonic scanning.

Model	SG-680
Table top	600 mm(W) * 1300 mm(L) 600 mm(W) * 1800 mm(L) with auxiliary table
Elevation	660 mm ~ 960 mm (without mattress)
Back / brake	0° ~ +70°
Tredelenbury	Head down/hip up 25° & head up/hip down 10°
Power supply	110V or 220V, 240V 50/60HZ. (can be followed by order requirement)



Standard accessories:

- | | |
|--|-------|
| 1. Pillow | 1 pc |
| 2. Waist side handle with clamps | 2 pcs |
| 3. Arm board | 2 pcs |
| 4. Leg crutches with clamps | 2 pcs |
| 5. Delivery collection plate with stainless steel filter | 1 pc |
| 6. Foot switch | 1 set |
| 7. Power cord | 1 pc |
| 8. Mattress | 1 set |
| 9. Hip pad (water proof) | 1 pc |
| 10. Auxiliary table with mattress | 1 pc |



Optional accessories:

- | | |
|--|-------|
| A. Head side hand holder with clamps | 1 set |
| B. Screen frame with clamps | 1 pc. |
| C. Heavy type leg crutch supporter with slipper and clamps | 2 pc. |

