

Gynecology-Delivery & operating SG-500 Series



SG-510:

Full function is controlled by Elec. Hydraulic system with memory elevation and memory back section up/down function moving range is adjustable.

SG-520:

Full function is controlled by Elec. Hydraulic system with memory elevation function; moving range is adjustable.

SG-530:

Full function is controlled by Elec. Hydraulic system without memory elevation function.

SG-540:

Elevation and back section are controlled by Elec. Hydraulic system. The seat section is immovable.



STURDY INDUSTRIAL CO., LTD.

MORE SAFETY - RELIABILITY-SATISFACTION

STURDY
APEX GROUP

Gynecology-Delivery & operating

SG-500 Series

Features:

- Fully automatic electric hydraulic control system.
- The memory device is built-in the elevation function. It can elevate while the memory up/down point automatically. (SG-520 only)
The memory device is built-in the elevation and back section function. It can be moved while the memories up/down point automatically. (SG-510 only)
- When the special position is needed, you can touch the foot control switch to adjust the back section up/down; table top up/down; and seat section up/down. (SG-510 & SG-520 & SG-530 only)
- The toilet paper roll placed under the seat. When the paper is properly set, you can get it continuously without moving the back or seat mattress.
- The surface of the table is covered with UV resistant plastic sheet, it is wear-resistant, pollution-proof and easy for maintenance.
- Water proof and 24V foot control switch.
- The back section can be moved to -15°, and the seat section can be moved +20°.(SG-540 is not available.) In emergency case it'll make the patients blood flow backward, to prevent the patient from suffering shock.
- It can be optioned the auxiliary table for delivery and operation use.
- All the accessories are made by stainless steel.



Specification:

TABLE TOP	55 cm (W) × 120 cm(L) 55 cm (W) × 175 cm(L)- WITH AUXILIARY
TABLE ELEVATION	60 cm ~85 cm (Without mattress) Maximum high/low can be followed by order requirement
BACK / BREAK	-15°~ +60°
WAIST / BREAK	0° TO +20° (SG-540 is not available)
Power Supply	110V or 220V/240V 50/60HZ (Can be followed by order requirement)

Standard accessories:

Screen Frame	1 PC.
Arm Board	1 PC.
Leg Crutches With Clamp	1 PR.
Forceps Cups	1 PC.
Foot Control Pedal	1 SET
Power Cord	1 PC.
Contamination Basin	1 PC.
Hip Pad (Water Proof)	1 PC.

Optional accessories:

Waist Side Handle With Clamp	2 PCS.
Delivery Collection Plate	1 PC.
Contamination Bucket, Tube & Funnel	1 SET
Auxiliary Table With Mattress	1 PC.
Heavy Type Leg Crutches With Clamps	2 PC.