



# Ultrasonic Cleaner

A professional ultrasonic cleaner is essential for removing the most stubborn residues from instruments.

Sturdy's ultrasonic cleaner ensure maximum efficiency and extreme ease of operation in everyday tasks and are available in a variety of models.

 Tanks Volume  
**0.6 - 22 Liters**



Reliable Sterilization Partner



Model No.	UC-400	UC-300	UC-150
Tank Size L x W x H (mm)	495 x 295 x 150	295 x 240 x 150	295 x 150 x 100
Tank Volume	22 Liters	10.5 Liters	4.5 Liters
Tank Material	#304 Stainless Steel	#304 Stainless Steel	#304 Stainless Steel
External Dimensions L x H x D (mm)	530 x 325 x 400	320 x 265 x 280	320 x 175 x 240
Housing Material	Aluminum Alloy (Powder Coating)	Aluminum Alloy (Powder Coating)	Aluminum Alloy (Powder Coating)
Lid Material	#304 Stainless Steel	#304 Stainless Steel	#304 Stainless Steel
Heating Power	400 W	400 W	200 W
Ultrasonic Power	400 W	300 W	150 W
Ultrasonic Frequency	40 kHz	46 kHz	46 kHz
Time Setting	1 - 99 mins	1 - 30 mins	1 - 30 mins
Temperature Setting	Room Temp. to 99 °C	Room Temp. to 60 °C	Room Temp. to 60 °C
Basket Material	#304 Stainless Steel	#304 Stainless Steel	#304 Stainless Steel
Power Supply	100 - 110 V / 220 - 240 V	100 - 110 V / 220 - 240 V	100 - 110 V / 220 - 240 V



Model No.	UC-80	UC-50
Tank Size L x W x H (mm)	150 x 135 x 100	168 x 89 x 55
Tank Volume	2 Liters	0.6 Liters
Tank Material	#304 Stainless Steel	#304 Stainless Steel
External Dimensions L x H x D (mm)	175 x 160 x 210	320 x 175 x 240
Housing Material	Aluminum Alloy (Powder Coating)	ABS Plastic
Lid Material	#304 Stainless Steel	Translucent PS Plastic
Heating Power	N/A	N/A
Ultrasonic Power	80 W	50 W
Ultrasonic Frequency	46 kHz	46 kHz
Time Setting	1 - 30 mins	Fixed
Temperature Setting	N/A	N/A
Basket Material	#304 Stainless Steel	ABS Plastic
Power Supply	100 - 110 V / 220 - 240 V	100 - 110 V / 220 - 240 V

**Sturdy**  
Welltell Group

No. 168, Sec. 1, Zhongxing Road, Wugu District,  
New Taipei City 248012, Taiwan  
T: +886 2 2981 5245  
F: +886 2 2987 8208 (Office)  
E: foreign@sturdy.com.tw  
[www.sturdy.com.tw](http://www.sturdy.com.tw)



Sturdy reserves the rights to improve the equipment shown without prior notice.

